

ISOQuant 1.6

Help pages

Jörg Kuharev, Pedro Navarro and Stefan Tenzer

September 29, 2016

NOTE:

Some contents of this document are not
aligned with the current software version.

An update is coming soon ...

Table of contents

1	Description	3
1.1	About ISOQuant	3
1.2	Citing ISOQuant	3
1.3	ISOQuant workflow	4
1.4	Known problems	6
2	Program requirements	7
3	Data requirements	8
3.1	Rawdata type	8
3.2	Database searches	8
3.3	Project design	8
3.3.1	Expression analysis	8
3.3.2	ISOQuant Project Designer	8
3.4	Peak detection/alignment/clustering	8
4	GUI and control elements	9
4.1	Main view	9
4.2	Project Finder	10
4.3	Context menu for PLGS projects in file system	11
4.4	Context menu for projects in database	12
4.5	Expression Analysis Selector	14
4.6	Project Designer	15
5	ISOQuant configuration	15
6	Configuration guide	17
6.1	EMRT cluster annotation	17
6.1.1	Peptide identification filter	17
6.1.2	Annotation mode	17
6.1.3	Annotation conflict filter	18
6.1.4	Homology / isoform and FDR filtering	18
6.2	Data preprocessing	18
6.3	EMRT table creation	18
6.4	Retention time alignment	18
6.5	EMRT clustering	19
6.6	Peak intensity normalization	19
6.7	Protein quantification	20
6.7.1	Peptide filtering	20
6.7.2	Protein quantification setting	20
6.8	Application behavior	20
6.8.1	User interface	20
6.8.2	Data source	21
6.8.3	Report	21
7	End-user license agreement	22
7.1	External components	22
7.2	ISOQuant license agreement	23
7.3	Disclaimer in plain language	23
8	About developers	24
	References	25

1 Description

ISOQuant is an academically developed, integrated bioinformatics pipeline for in-depth evaluation and statistical data analysis of data-independent acquisition (MS^E and IMS-MS^E) based label-free quantitative proteomics that improves data reliability and quality by application of well-established and novel analysis methods.

1.1 About ISOQuant

One of the main bottlenecks in the evaluation of label-free quantitative proteomics experiments is the often cumbersome data export for in-depth data evaluation and analysis. Data-independent, alternate scanning LC-MS (MS^E /HDMS E /UDMS E) peptide fragmentation data can currently only be processed by Waters PLGS software (and by recently introduced Progenesis QI for Proteomics (Nonlinear Dynamics / Waters)).

PLGS performs absolute quantification only on a run-to-run level, it does not afford absolute quantification of protein isoforms and label-free relative quantification of peptides and proteins based on clustered accurate mass-retention time pairs on a complete experiment basis.

The bioinformatics pipeline ISOQuant directly accesses xml files from the PLGS root folder and browses for relevant data from a label-free Expression Analysis project (quantification analyses, sample descriptions, search results) for fully automated import into a MySQL database. EMRTs are subjected to multidimensional LOWESS-based intensity normalization and annotated by matching exact masses and aligned retention times of detected features with highest scoring peptide identification data from associated workflows. Based on the annotated cluster table, ISOQuant calculates absolute in-sample amounts with an integrated protein isoform quantification method, utilizing average intensities of proteotypic peptides for the partitioning of non-unique peptide intensities between protein isoforms. All data is stored in a local MySQL based database that can be queried directly by experienced users.

1.2 Citing ISOQuant

ISOQuant has been developed since 2009. We introduced the basic principles of ISOQuant analysis to the community as part of a Nature Methods articles 2014 (Distler et al. 2014).

```
Distler, U., Kuharev, J., Navarro, P., Levin, Y., Schild, H., & Tenzer, S.
(2014). Drift time-specific collision energies enable deep-coverage
data-independent acquisition proteomics.
Nature Methods, 11(2), -167170. http://doi.org/10.1038/nmeth.2767
```

Please cite the mentioned publication, when using ISOQuant to produce publication data or referencing to ISOQuant in other context. Use the following BibTeX code to import into the reference manager of your choice:

```
@article{distler_drift_2014,
  title = {Drift time-specific collision energies enable
    deep-coverage data-independent acquisition proteomics},
  author = {Distler, Ute and Kuharev, Jörg and Navarro, Pedro
    and Levin, Yishai and Schild, Hansjörg and Tenzer, Stefan},
  journal = {Nature Methods},
  volume = {11},
  issn = {1548-7091},
  url = {http://www.nature.com/nmeth/journal/v11/n2/full/nmeth.2767.html},
  doi = {10.1038/nmeth.2767},
  month = feb,
  year = {2014},
  pages = {167--170}
}
```

1.3 ISOQuant workflow

The data analysis workflow (see fig. 1) consists of raw data preprocessing using vendor software PLGS and the downstream analysis using ISOQuant. In our data analysis workflow, PLGS is used for the initial signal processing as well as for peptide and protein identification. Before automatically importing PLGS results into a relational database (MySQL), ISOQuant allows to change the structure of underlying PLGS project and in this way to redesign the label-free experiment. The ISOQuant data analysis workflow is composed of multiple dedicated algorithms. At different stages of analysis, data is filtered on peptide and protein level based on user defined criteria (identification score and type, sequence length, replication rate, FDR threshold, etc.) to ensure a constant level of high data quality for all runs in the project. The retention time alignment procedure corrects non-linear retention time distortions between LC-MS runs of the experiment (Podwojski et al. 2009). To group corresponding features from different runs of the experiment, exact mass and retention time pairs (EMRT) extended by ion mobility values are evaluated using the density based clustering algorithm DBSCAN (Ester et al. 1996). Resulting feature clusters are annotated by evaluation of consent peptide identifications and identification probabilities. The feature cluster annotation approach transfers peptide identifications between runs and reduces missing values increasing the reproducibility of data analysis. Resolving ambiguous peptides-in-proteins networks, the protein homology filtering algorithm reduces side effects of the protein inference problem (Nesvizhskii and Aebersold 2005). The multi-dimensional feature intensity normalization algorithm reduces the effects of technical variance between independently acquired LC-MS runs of the experiment and increases the reproducibility of quantification. Finally, ISOQuant uses a derivative of the Top3 method (Silva et al. 2006) for the absolute protein quantification and exports analysis results to multiple output formats for subsequent interpretation. Alternatively to ISOQuant, the commercial software Progenesis QI for proteomics (Nonlinear Dynamics / Waters) or the freely available R-package synapter (Bond et al. 2013) could also be used to analyze MS^E data. In a recent study, we tested Progenesis QIP, synapter and ISOQuant for the performance of protein identification and the label-free quantification based on the analysis a complex metaproteome sample set (Kuharev et al. 2015).




Figure 1: The workflow of ISOQuant analysis

1.4 Known problems

- On some **Windows Vista** or **Window 7** machines ISOQuant can not write its configuration file. In this case you have to execute ISOQuant with administrative privileges or correct file system permissions for ISOQuant installation folder. This is not an ISOQuant issue, sometimes Windows messes up file system permissions by using different and inconsistent user privileges at different time points.
- Analysis of high complexity datasets may take a while.
- Importing (and analyzing) large projects or runs of high complexity may cause out of memory errors, make sure your PC has enough memory and assign more Heap-Space to Java Virtual Machine for running ISOQuant Application.
- ISOQuant may fail to import and process data if some PLGS project files are broken.
- Running MySQL on Mac OSX machines significantly decreases the performance of ISOQuant. This is a known problem of MySQL not of ISOQuant. Use Windows or Linux machines and/or install MariaDB instead of MySQL for better performance.

2 Program requirements

ISOQuant will only work properly if the system for running ISOQuant meets following requirements.

- Operating System: Windows, Mac OS X or Linux
- PLGS root folder with projects containing processed MS^E/HDMS^E/UDMS^E data is accessible (tested PLGS versions: 2.3/2.4/2.5/3.0)
- at least 3GB RAM
- Java Runtime Environment version 1.6.0 (or newer) is installed and works properly
- MySQL Server 5.1 (or newer) is installed and running on local machine or network. (tested MySQL versions: 5.1 - 5.5)
- MySQL configuration file options for heap and temporary tables have large values as shown in following listing. In some cases it may be useful also to increase the size of MySQL thread stack.

```
max_heap_table_size = 2048M
tmp_table_size = 2048M
thread_stack = 256K
```

Depending on your operating system and MySQL-Version the configuration file is named either my.ini or my.cnf and its location may vary.

Following listing shows an example of MySQL configuration section [mysqld] working for us on MacOSX 10.6.8 Snow Leopard running MySQL Server from XAMPP 1.7.3, the configuration file is located in /Applications/XAMPP/xamppfiles/etc/my.cnf

```
[mysqld]
port = 3306
socket = /Applications/XAMPP/xamppfiles/var/mysql/mysql.sock
skip-locking
key_buffer = 128M
max_allowed_packet = 16M
table_cache = 128
sort_buffer_size = 32M
read_buffer_size = 8M
read_rnd_buffer_size = 8M
net_buffer_length = 64K
thread_stack = 256K
myisam_sort_buffer_size = 32M
tmpdir = /Applications/XAMPP/xamppfiles/temp/
max_heap_table_size = 2048M
tmp_table_size = 2048M
sync_frm = 0
skip-sync_frm=OFF
```

Do not forget to restart MySQL after editing its configuration!

Expert note:

If you get some “out of memory” errors while running ISOQuant please make sure you start the application by giving Java Virtual Machine a chance to have enough memory space by command line options, e.g.

```
java -Xms256m -Xmx2G -jar ISOQuant.jar
```

this command will assign up to 2 GBs (parameter **-Xmx2G**) RAM to the virtual machine and run the ISOQuant application. For some very complex datasets, it may be useful to increase the **-Xmx** value e.g. **-Xmx48G** to allow the virtual machine to access 48 GBs of RAM.

3 Data requirements

3.1 Rawdata type

ISOQuant has been developed for Waters QTOF LC-MSE and Waters Synapt G2/G2-S LC-MS^E/HDMS^E/UDMS^E instrument data. At this time, only 1D-UPLC data is fully supported, 2D-UPLC support will be included in later releases.

3.2 Database searches

At the moment, ISOQuant can only process Ion Accounting workflows (MS^E/HDMS^E/UDMS^E-data). Classical DDA-Type experiments are not yet supported.

3.3 Project design

There are two different ways to use ISOQuant either as an extension to PLGS Expression Analysis or completely replacing it.

3.3.1 Expression analysis

You can use your experiment design given by running PLGS Expression analysis. As a prerequisite for this approach, a complete expression analysis of multiple samples and replicates is required. In PLGS, please select autonormalization of samples for generating EMRT and protein tables. Both EMRT and Protein tables have to be created during the expression analysis. Each expression analysis within a PLGS project can be selected during processing in ISOQuant. ISOQuant will create separate databases for storing the data of every single processed expression analysis.

3.3.2 ISOQuant Project Designer

As an alternative to the PLGS Expression analysis, you can use the simple and efficient built-in Project Designer described in section 4.6.

3.4 Peak detection/alignment/clustering

ISOQuant is based on peak detection, alignment and clustering of data performed by PLGS. We are aware of some peak splitting/alignment/clustering issues in PLGS. Therefore, we have spent a lot of time to develop own methods for these tasks. You can either keep EMRT alignment/clustering results or let ISOQuant do the complete analysis. For details see future publications.

4 GUI and control elements

4.1 Main view

Figure 2 shows the main view of ISOQuant. User interaction is applied by following control elements:

1. List of projects found in PLGS root folder
2. List of projects from ISOQuant database
3. Button choose PLGS root folder
4. Button choose database
5. Button restore project from file
6. Button find projects
7. Button edit configuration
8. Button show help window
9. Button shutdown application




Figure 2: the main view of ISOQuant

4.2 Project Finder

The Project Finder window as shown in figure 3 makes possible to search projects lists for projects by substrings of project titles and regular expressions matching to them. In case your search string matches to one or multiple projects the Project Finder will mark these projects by selecting them in both file system and database projects lists. The Project Finder window can be accessed by clicking the button find projects from the tool bar on the main application window.




Figure 3: Project Finder window

4.3 Context menu for PLGS projects in file system

Advanced options for each project from PLGS root folder are available from a context menu like shown in figure 4:

1. find in database
finds selected projects in the list of projects from database by comparing their titles and select them if such projects exist.
2. about project
shows additional information about selected projects.
3. import and process
allows to select one of predefined processing queues and starts processing selected projects using selected processing queue.




Figure 4: Context menu for PLGS projects in file system

4.4 Context menu for projects in database

Advanced options for each project from database are available from a context menu like shown in figure 5

1. find in file system
finds selected projects in the list of projects from file system by comparing their titles and select them if such projects exist.
2. show info
shows additional information about selected projects.
3. rename project
rename selected projects.
4. reprocess
reprocess a project starting from user selected processing stage. All needed subsequent processing steps are automatically applied.
5. create report
generate on of implemented report types.
6. export to file
export selected projects from database to (backup) files which can be imported by other ISOQuant instances.
7. remove from database
removes selected projects from database.

Note: We continue to develop and improve ISOQuant. The context menus could look different in different software releases. The elements in context menus are subject of change, their number and order may vary.




Figure 5: Context menu for already processed projects

4.5 Expression Analysis Selector

In some cases a single PLGS project contains multiple defined Expression Analyses. Some processing queues work with project structures provided by PLGS Expression Analyses. These queues start with the selection of contained expression analyses for each selected project. The selection of expression analyses is done within the Expression Analysis Selector window as shown in figure 6 by activating the checkboxes from the column include for each Expression Analysis to be processed. The Expression Analysis Selector shows each previously selected project in its own tab pane. ISOQuant generates a separate database for each selected Expression Analysis.




Figure 6: Expression Analysis Selector

4.6 Project Designer

ISOQuant allows to create user defined project structures and then to process these newly structured project data. User defined project structures are created using the Project Designer as shown in figure 7. The Project Designer window shows the PLGS project structure on the left and the user defined structure on the right. A new project structure is created by drag and drop based moving of workflows, samples or groups between left and right structure trees. Additionally to drag and drop actions, right click context menus are available on the right side of Project Designer enabling editing and removing of selected structure elements. On the top of window you can switch between Project Designer panes to restructure each previously selected PLGS project. Processing of designed project can be initiated by clicking the button Ok on the bottom of window or can be aborted by clicking the Cancel button. While processing ISOQuant generates a separate database for each designed project.




Figure 7: Project Designer

5 ISOQuant configuration

ISOQuant stores parameters for program behavior and data processing algorithms in a single configuration file named **isoquant.ini**. This configuration file is located in the folder you have installed ISOQuant to. For resetting parameters to default values just close ISOQuant then delete or rename the configuration file (or single parameter lines) and start ISOQuant again. If no configuration file can be found on application start a new one will be created using default parameter values.

Do not change the configuration file unless you know what you do!

ISOQuant configuration can be edited from ISOQuant Configuration Editor accessible from graphical user interface. Configuration Editor allows to edit parameters and also export/import configuration settings to/from files.

Two configuration files are provided with ISOQuant installation packages:

- **isoquant_high_confidence.ini** example configuration file for high confidence quantitative analyses
- **isoquant_maxID.ini** example configuration file for discovery proteomics experiments

These files can be imported into ISOQuant from Configuration Editor or manually copied to **isoquant.ini** file.

6 Configuration guide

This chapter lists and describes the main set of available parameters. The number of parameters, their names and the behavior of application caused by parameters are subjects of change because we actively work on improving ISOQuant and underlying methods. Thus the following list of parameters may be incomplete and/or up to date.

6.1 EMRT cluster annotation

6.1.1 Peptide identification filter

User can define minimum reliability criteria for a peptide to be used for ISOQuant processing. *Note:* Only peptides passing this filter will be used for all further analysis steps!

6.1.1.1 Peptide types

Peptides identified in first pass of PLGS database search (type: **PEP_FRAG_1**) are generally accepted. Setting one of following parameters to TRUE will configure ISOQuant to accept additional peptide types.

- **process.identification.peptide.acceptType.IN_SOURCE**=false
- **process.identification.peptide.acceptType.MISSING_CLEAVAGE**=false
- **process.identification.peptide.acceptType.NEUTRAL_LOSS_H2O**=false
- **process.identification.peptide.acceptType.NEUTRAL_LOSS_NH3**=false
- **process.identification.peptide.acceptType.PEP_FRAG_2**=false
- **process.identification.peptide.acceptType.PTM**=false
- **process.identification.peptide.acceptType.VAR_MOD**=false

Following additional filtering criteria are used as thresholds to select peptides used for further processing steps.

- **process.identification.peptide.minReplicationRate**=2.0
minimum acceptable peptide replication rate based on absolute number of runs in which every peptide (as sequence-modifier tuple) was identified.
- **process.identification.peptide.minScore**=1.0
minimum acceptable PLGS peptide identification score.
- **process.identification.peptide.minOverallMaxScore**=1.0
minimum acceptable value of highest PLGS identification score of a peptide reached in any run of a project. For a peptide detected in multiple runs its maximum reached score has to hit this score to be accepted for annotation. Increasing score will reduce the number of peptides used for EMRT cluster annotation and not necessarily the overall protein quantification quality. Recommended values are between 0.0 and 5.0
- **process.identification.peptide.minSequenceLength**=6
minimum acceptable peptide sequence length. Recommended value is 6 or more.

6.1.2 Annotation mode

- **process.annotation.useSharedPeptides**=all
 - **all** this is the normal case.
 - **unique** only unique peptides are used for further processing, this option removes all shared peptides from peptides-in-proteins relation instead of protein homology filtering solving the problem of protein inference in a very radical way.
 - **razor** only razor and unique peptides are used for further processing, this option removes all shared peptides from peptides-in-proteins relation after protein homology filtering. Razor and unique peptides are highly reliable for protein quantification because their intensity can be directly assigned to a protein.

6.1.3 Annotation conflict filter

There are cases when multiple peptide identifications map to a single EMRT cluster.

- **process.annotation.peptide.sequence.maxCountPerEMRTCluster=1**
acceptable number of different peptide sequences (remaining after filtering peptides) allowed to annotate a single cluster. The annotation process will skip ambiguous clusters if this value is set to **1**. For bigger values, the annotation conflicts are resolved by annotating clusters with the peptide having the highest sum of PLGS identification scores in this cluster.

6.1.4 Homology / isoform and FDR filtering

- **process.annotation.protein.resolveHomology=true**
Should proteins be filtered for (peptide sequence based) homology/isoform. Only one of detected homologue proteins will be reported.
- **process.annotation.peptide.maxFDR=0.01**
maximum accepted false discovery rate level for peptides. (value 0.01 means maximum 1% FDR)

6.2 Data preprocessing

- **process.peptide.deplete.PEP_FRAG_2=false**
should PEP_FRAG_2 peptides be completely removed from database.
- **process.peptide.deplete.CURATED_0=false**
should CURATED=0 peptides be completely removed from database. If **true**, low-quality peptide IDs are removed.

6.3 EMRT table creation

- **process.emrt.minIntensity=1000**
peaks having intensities below this limit are assumed to be noise and will not appear in EMRT table
- **process.emrt.minMass=500**
peaks having masses below this limit are assumed to be noise and will not appear in EMRT table

6.4 Retention time alignment

- **process.emrt.rt.alignment.match.maxDeltaMass.ppm=10.0**
maximum accepted mass difference between two signals to be assumed as matching for retention time alignment.
- **process.emrt.rt.alignment.match.maxDeltaDriftTime=2.0**
maximum accepted drift time (ion mobility) difference between two signals to be assumed as matching for retention time alignment. This value is ignored for non-ion-mobility projects. Large value, e.g. 200 will disable the effect of ion mobility on the time alignment.
- **process.emrt.rt.alignment.minIntensity=1000**
only peaks with intensity over this threshold value are considered for the retention time alignment procedure
- **process.emrt.rt.alignment.minMass=800.0**
only peaks with mass over this threshold value are considered for the retention time alignment procedure
- **process.emrt.rt.alignment.normalizeReferenceTime=false**
if true, resulting reference times are adjusted to median distortions at every time point.
- **process.emrt.rt.alignment.maxProcesses=4**
the maximum number of concurrent retention time alignment processes. We recommend values between 1 and the number of available CPU cores. Default value is set to $\frac{1}{2}$ of the number of CPU cores
- **process.emrt.rt.alignment.maxProcessForkingDepth=4**
the maximum multithreading depth for each retention time alignment process

6.5 EMRT clustering

- **process.emrt.clustering.distance.unit.mass.ppm=6.0**
minimum mass based distance between clusters. This is an instrument dependent parameter, e.g. 6 ppm is a good value for Waters Synapt G2/G2-S and 10-12 ppm for Waters Q-TOF Premier and Synapt G1
- **process.emrt.clustering.distance.unit.time.min=0.2**
minimum retention time based distance between clusters. This is a LC gradient length and peak width dependent parameter, good values are observed to be between 0.06 and 0.2, we recommend to try 0.08, 0.12, 0.16, 0.2; please report which values would work for your setup at which gradient length
- **process.emrt.clustering.distance.unit.drift.bin=2.0**
minimum drift time based distance between clusters. This is an instrument and also project setup dependent parameter, e.g. for pure IMS projects containing G2 or G2S data, we recommend a value of 2.0. This value is ignored for non-ion-mobility projects. Large value, e.g. 200 will disable the effect of ion mobility on the EMRT clustering.
- **process.emrt.clustering.dbscan.minNeighborCount=2**
the minimum cluster size (except of noise) and also the minimum required number of peaks inside the reachability radius for cluster expansion. This is a DBSCAN specific parameter and should be increased for big projects. The value of this parameter also depends on used clustering distance units.
- **process.emrt.clustering.maxProcesses=8**
the maximum number of concurrent clustering processes. For best performance is reached by setting this value to the number of available CPU cores. Default value is set by the estimated number of available CPU cores.

6.6 Peak intensity normalization

- **process.normalization.minIntensity=3000**
systematic errors of peptides with intensities below this limit are ignored during normalization process.
- **process.normalization.lowess.bandwidth=0.3**
bandwidth parameter for non-linear regression method (LOWESS) used for exploring systematic errors during normalization process. Recommended values are between 0.3 and 0.6
- **process.normalization.orderSequence=XPIR**
The processing order sequence of dynamic multi-dimensional normalization. The processing order sequence is defined as a word build from following characters: **X, P, G, S, W, I, R, M, E**. The occurrence of a letter either defines the next dimension for EMRT normalization or changes the normalization mode or discards previously calculated values:
 - X** reset emrt intensities to original values
 - P** activate IN-PROJECT normalization mode, average intensity of all emrts in a cluster is used as the reference.
 - G** activate IN-GROUP normalization mode, average intensity of all emrts from a group of samples within a cluster is used as the reference. T.m. each run uses reference values from its sample group.
 - S** activate IN-SAMPLE normalization mode, average intensity of all emrts from a sample within a cluster is used as the reference. T.m. each run uses reference values from its sample.
 - W** activate Workflow/Run-Value based normalization mode, the run to be the normalization reference is automatically set by choosing the run with the highest number of emrts.
 - I** normalize emrt intensities using log-intensity dimension
 - R** normalize emrt intensities using retention time dimension
 - M** normalize emrt intensities using mass dimension
 - E** equalize emrt intensities by adjusting sums of intensities for each run

The order sequence is processed from left to right, e.g. the recommended order sequence **XPIR** stands for clean in-project normalization using intensity domain followed by normalization using retention time domain.

6.7 Protein quantification

6.7.1 Peptide filtering

Peptides for protein quantification may be filtered by their type and minimum reached score of a peptide per EMRT cluster. **PEP_FRAG_1** peptides are always accepted. User may decide to accept additional peptide types for quantification. Allowing additional types may result in higher number of quantified proteins but also may affect the quality of quantification. *Note:* This peptide filtering step can not recover peptides not passed the peptide identification filter.

- **process.quantification.peptide.acceptType.IN_SOURCE**=false
- **process.quantification.peptide.acceptType.MISSING_CLEAVAGE**=false
- **process.quantification.peptide.acceptType.NEUTRAL_LOSS_H2O**=false
- **process.quantification.peptide.acceptType.NEUTRAL_LOSS_NH3**=false
- **process.quantification.peptide.acceptType.PEP_FRAG_2**=false
- **process.quantification.peptide.acceptType.PTM**=false
- **process.quantification.peptide.acceptType.VAR_MOD**=false
- **process.quantification.peptide.minMaxScorePerCluster**=5.0

6.7.2 Protein quantification setting

- **process.quantification.absolute.standard.entry**=ENO1_YEAST
entry of protein used as quantification standard
- **process.quantification.absolute.standard.fmol**=50.0
amount of quantification standard protein
- **process.quantification.absolute.standard.used**=true
is a quantification standard protein used at all?
- **process.quantification.topx.degree**=3
maximum number of peptides for quantifying single proteins
- **process.quantification.maxProteinFDR**=0.01
maximum accepted false discovery rate level for reported proteins. (value 0.01 means 1% FDR level)
- **process.quantification.minPeptidesPerProtein**=1
a protein is reported only if it can be quantified by using as minimum this number of peptides

6.8 Application behavior

6.8.1 User interface

- **setup.ui.captureConsoleMessages**=true
show Java console messages inside ISOQuant message panel
- **setup.ui.location.left**=560
ISOQuant window location, pixels from left
- **setup.ui.location.top**=360
ISOQuant window location, pixels from top
- **setup.ui.size.height**=480
ISOQuant window height
- **setup.ui.size.width**=800
ISOQuant window width
- **setup.ui.promptForExit**=true
ask user on closing window
- **setup.ui.iconScaleFactor**=1.0
scale original icon sizes by this factor, may be useful on unusually small or large screens.

6.8.2 Data source

- **setup.db.autoLoad=false**
should application connect database on start
- **setup.db.host=localhost**
the MySQL database host name or ip address. Default tcp port number for MySQL servers is 3306. If your MySQL server is running using an other tcp port number, its host name has to be expanded by adding ':' and the correct port number, e.g. localhost:3307 for MySQL server running on local machine and listening at tcp port 3307
- **setup.db.user=root**
MySQL user name
- **setup.db.pass=**
MySQL users password
- **setup.plgs.root.showEACount=true**
should number of PLGS expression analyses be determined and shown
- **setup.plgs.root.showFSSize=false**
should file system size of a projects folder be determined and shown
- **setup.plgs.root.dir=/Volumes/RAID0/PLGS2.5/root**
path of last selected PLGS root folder
- **setup.plgs.root.autoLoad=false**
should application read last used root folder on start

6.8.3 Report

- **setup.report.dir=/Volumes/RAID0/reports**
path of last selected report output folder.
- **setup.report.csv.columnSeparator=','**
column separator string (enclosed in ' or "), usually either ',' or ';'.
- **setup.report.csv.decimalPoint=''**
decimal point string (enclosed in ' or "), usually either '' or ' '.
- **setup.report.csv.textQuote='"**
string for quoting text blocks, usually
- **setup.report.mzidentml.DBNCBITaxID=**
- **setup.report.mzidentml.DBOrganismScientificName=**
- **setup.report.mzidentml.DBversion=**
- **setup.report.mzidentml.researcherFirstName=John**
- **setup.report.mzidentml.researcherLastName=Doe**
- **setup.report.mzidentml.researcherOrganization=Uni-Mainz**
- **setup.report.xls.showAbsQuantFMOLUG=true**
create an extra sheet for absolute protein quantity in femtomoles per microgram.
- **setup.report.xls.showAbsQuantFMOL=true**
create an extra sheet for absolute protein quantity in femtomoles.
- **setup.report.xls.showAbsQuantNG=true**
create an extra sheet for absolute protein quantity in nanograms.
- **setup.report.xls.showAbsQuantPPM=true**
create an extra sheet for absolute protein quantity in parts per million.
- **setup.report.xls.showAllProteins=false**
create an extra sheet for some PLGS based proteins details.
- **setup.report.xls.showRTAlignment=false**
create an extra sheet for retention alignment results.

Generated excel sheets are limited to maximum 65536 rows when using old XLS (Excel 97/2000/2003) format due to its technical limitations. In case of doubt, please use XLSX format for creating ISOQuant reports.

7 End-user license agreement

7.1 External components

ISOQuant relies on external components (re)distributed under different conditions. By using ISOQuant you need to agree to terms and conditions of third party libraries and software included in ISOQuant or needed for running ISOQuant.

ISOQuant uses following external Java libraries:

Library	Version	License Type	Purpose
JDOM	1.1.3	BSD License	XML handling
Tagsoup	1.2.1	Apache v2.0	XML parsing
MySQL Connector/J	5.1.13	GPLv2 with FOSS Exception	database communication
JSiX	1.0	BSD License	Java extensions
POI	3.8	Apache v2.0	spreadsheet file creation
DOM4J	1.6.1	BSD License	POI dependency
StAX	1.0.1	Apache v2.0	POI dependency

Binary versions of these libraries are repackaged into and redistributed with ISOQuant software package according to their license conditions. Please find original licenses as part of ISOQuant package.

Furthermore, ISOQuant relies on external environmental software being not a part or a component of ISOQuant but needed to run it:

- Operating System
- Java Virtual Machine
- MySQL 5.1 compatible database engine (e.g. MySQL Server: <http://www.mysql.com/> or MariaDB: <http://mariadb.org/>)
- Waters ProteinLynx Global Server

Please pay attention to terms and conditions arising from any software usage in any way related to ISOQuant.

7.2 ISOQuant license agreement

ISOQuant binaries and source code are available under conditions of BSD license as follows:

```
ISOQuant - integrated solution for LC-MS based label-free protein quantification

Copyright (c) 2009 - 2013, JOERG KUHAREV and STEFAN TENZER
All rights reserved.

Redistribution and use in source and binary forms, with or without
modification, are permitted provided that the following conditions are met:
1. Redistributions of source code must retain the above copyright
   notice, this list of conditions and the following disclaimer.
2. Redistributions in binary form must reproduce the above copyright
   notice, this list of conditions and the following disclaimer in the
   documentation and/or other materials provided with the distribution.
3. All advertising materials mentioning features or use of this software
   must display the following acknowledgement:
   This product includes software developed by JOERG KUHAREV and STEFAN TENZER.
4. Neither the name "ISOQuant" nor the
   names of its contributors may be used to endorse or promote products
   derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY JOERG KUHAREV and STEFAN TENZER 'AS IS' AND ANY
EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED
WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE
DISCLAIMED. IN NO EVENT SHALL JOERG KUHAREV BE LIABLE FOR ANY
DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES
(INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES;
LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND
ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT
(INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS
SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.
```

7.3 Disclaimer in plain language

- ISOQuant is an academically developed freely available software implemented using only freely available technologies
- By using ISOQuant you automatically agree to terms and conditions of third party libraries and software included in ISOQuant or needed for running ISOQuant
- We provide ISOQuant as is, without any kind of warranty
- You are using ISOQuant at your own risk, we are not responsible for any kind of damage, errors or data loss produced in any way related to ISOQuant.

8 About developers

ISOQuant is developed by

- Jörg Kuharev (kuharev@uni-mainz.de) and
- Pedro Navarro (pnavarro@uni-mainz.de) and
- Stefan Tenzer (tenzer@uni-mainz.de)

UNIVERSITÄTSMEDIZIN
der Johannes Gutenberg-Universität
Institut für Immunologie
Core Facility für Massenspektrometrie
Germany, Mainz, September 29, 2016

References

- Bond, Nicholas J, Pavel V Shliha, Kathryn S Lilley, and Laurent Gatto. 2013. “Improving Qualitative and Quantitative Performance for MS E-Based Label-Free Proteomics.” *Journal of Proteome Research*.
- Distler, Ute, Jörg Kuharev, Pedro Navarro, Yishai Levin, Hansjörg Schild, and Stefan Tenzer. 2014. “Drift Time-Specific Collision Energies Enable Deep-Coverage Data-Independent Acquisition Proteomics.” *Nature Methods* 11 (2): 167–70. doi:10.1038/nmeth.2767.
- Ester, Martin, Hans P. Kriegel, Jörg Sander, and Xiaowei Xu. 1996. “A Density-Based Algorithm for Discovering Clusters in Large Spatial Databases with Noise.” In *Second International Conference on Knowledge Discovery and Data Mining*, edited by Evangelos Simoudis, Jiawei Han, and Usama Fayyad, 226–31. Portland, Oregon: AAAI Press.
- Kuharev, Jörg, Pedro Navarro, Ute Distler, Olaf Jahn, and Stefan Tenzer. 2015. “In-Depth Evaluation of Software Tools for Data-Independent Acquisition Based Label-Free Quantification.” *PROTEOMICS*. doi:10.1002/pmic.201400396.
- Nesvizhskii, Alexey I., and Ruedi Aebersold. 2005. “Interpretation of Shotgun Proteomic Data: The Protein Inference Problem.” *Molecular & Cellular Proteomics: MCP* 4 (10): 1419–40. doi:10.1074/mcp.R500012-MCP200.
- Podwojski, Katharina, Arno Fritsch, Daniel C. Chamrad, Wolfgang Paul, Barbara Sitek, Kai Stühler, Petra Mutzel, et al. 2009. “Retention Time Alignment Algorithms for LC/MS Data Must Consider Non-Linear Shifts.” *Bioinformatics* 25 (6): 758–64. doi:10.1093/bioinformatics/btp052.
- Silva, Jeffrey C, Marc V Gorenstein, Guo-Zhong Li, Johannes P C Vissers, and Scott J Geromanos. 2006. “Absolute Quantification of Proteins by LCMSE: A Virtue of Parallel MS Acquisition.” *Molecular & Cellular Proteomics: MCP* 5 (1): 144–56. doi:10.1074/mcp.M500230-MCP200.